

IMMUNOSTAINING OF *Drosophila* NEUROBLASTS

- 1- Dissect the brains (without imaginal discs) of third instar larvae in 0.7% NaCl.
- 2- Fix the brains in 4% paraformaldehyde (or 3.7% formaldehyde) in cytoskeleton buffer (or PBS) for 30 min.
- 3- Incubate 3 min in 45% acetic acid (for γ -tubulin immunostaining immerse brains in methanol for 2 min before fixation in acetic acid).
- 4- Strongly squash the brains in 60% acetic acid under a coverslip and freeze in liquid nitrogen.
- 5- Carefully remove the coverslips and immediately immerse the slides for 10 min in ethanol at -20°C (inside a black box in a dry ice container).
- 6- Incubate for 10 min in PBS + 0.1% Triton (black box).
- 7- Wash 2x 5 min in PBS (black box).
- 8- Incubate for 45 min in 1% BSA in PBS for blocking (black box).
- 9- Incubate with primary antibodies in 1% BSA + 0.1% Triton in PBS for 1 hour.
- 10- Wash 3x 5 min in PBS + 0.1% Triton (black box).
- 11- Incubate with appropriate secondary antibodies in 1% BSA + 0.1% Triton in PBS for 1 hour.
- 12- Wash 3x 5 min in PBS + 0.1% Triton (black box).
- 13- Wash 1x 5 min in PBS (black box).
- 14- Mount the brains with DAPI diluted 1:1000 in Vecta-Shield and seal the coverslip.